


MASTERPULS® MP200



STORZ MEDICAL



MASTERPULS® MP200

The **MOST RELIABLE** shock wave generation

When launched in October 2003 as the first-generation MASTERPULS® system, the MASTERPULS® MP100 set new standards in radial shock wave therapy.

Successful, reliable, safe, versatile or »simply better« – these are the most frequently cited attributes of the MASTERPULS® series. Ease of use and simple handpiece servicing by the user at reduced costs are the key benefits of the MASTERPULS® MP50, MP100 and MP200 models.

The terms that best describe the MASTERPULS® product family today are innovation, safety and reliability! In the field of pain therapy and rehabilitation MASTERPULS® shock wave systems are a must-have for any modern, economically successful medical practice.

Indications	Success rate
Calcified tendinitis of the shoulder	> 83%
Radial and ulnar humeral epicondylitis	> 80%
Patellar tendinitis/Achillodynia	> 88%
Plantar fasciitis	> 80%
Myofascial trigger points	> 90%
Acupuncture shock wave therapy	up to 90%
Other insertional disorders	
Muscle and connective tissue activation with V-ACTOR®	



Key features at a glance

Combinable handpieces

- ▶ Radial shock wave therapy with various shock transmitters
- ▶ Vibration therapy (V-ACTOR®)

Pulse frequency/pressure

- ▶ Radial shock wave therapy: 1 – 21 Hz/1 – 5 bar
- ▶ Vibration therapy (V-ACTOR®): 1 – 35 Hz/1 – 5 bar

Oscillating »D-ACTOR®« technology

- ▶ for better myofascial trigger point therapy

»Deep Impact®« shock transmitter

- ▶ for the treatment of deep pain regions

»CERAmA-x™« shock transmitter

- ▶ elastic shock transmitter for shock wave application without gel

»V-ACTOR®« handpiece

- ▶ for muscle and connective tissue activation/smoothing