





# Shockwave

## CANADA<sup>INC</sup>

	MASTERPULS® ONE	MASTERPULS® R-SW50	MASTERPULS® ultra+	MASTERPULS® 200
Size	 23.9cm/28.9cm/31cm	 26.2cm/18.6cm/24.2cm	 42.6cm/14.4cm/34cm	 46.6cm/18.7 cm/45.4cm
Weight	21 lbs	15.4 lbs	21 lbs	50 lbs
Power Level	2.7 Bar	4 Bar	5 Bar	5 Bar
Frequency	6-18 Hz	1-18 Hz	1-21 Hz	1-21 Hz
Sparrow Handpiece	✓			
Falcon Handpiece		✓	✓	✓
C-15,D20-S	✓			
C-15, DI15, D20-T			✓	✓
Complete set of 9 transmitter heads		Compatible	Compatible	Compatible
Storz Tablet		Compatible (Optional)		Included
Vacu-Actor cupping technology and materials				Included
V-Actor Handpiece		Compatible	Compatible	Compatible
Set of Fascia Transmitter Heads		Compatible	Compatible	Compatible
Set of Spinal Transmitter Heads		Compatible	Compatible	Compatible



# MP ONE

SIZE (WIDTH/HEIGHT/DEPTH)	KEY FEATURES	COMPATIBILITY
<p>23.9cm/28.9cm/31cm</p> <p>21lbs</p>	<p>Easy to Carry Handle</p> <p>Touch Screen</p> <p>Sparrow Handpiece</p> <p>6 Energy Levels</p> <p>6-18 HZ</p> <p>Max 2.7 Bar</p>	<p>C15, DI15, D20-S Transmitter Heads</p>



# MP R-SW50

SIZE (WIDTH/HEIGHT/DEPTH)	KEY FEATURES	COMPATIBILITY
<p>26.2cm/18.6cm/24.2cm</p> <p>15.4lbs</p>	<p>Patented Handpiece Control Technology</p> <p>Falcon Handpiece</p> <p>1-18 HZ</p> <p>Max 4.0 Bar</p>	<p>Complete set of 9 transmitter heads</p> <p>Storz Tablet</p> <p>V-Actor Handpiece</p> <p>Fascia Transmitter Heads</p> <p>Spinal Transmitter Heads</p> <p>Storz Cart</p>

Please note that optional accessories are highlighted in red text. Everything in black text is included with the device.



# MP ultra+

SIZE (WIDTH/HEIGHT/DEPTH)	KEY FEATURES	COMPATIBILITY
<p>42.6cm/14.4cm/34cm</p> <p>21lbs</p>	<p>Falcon Handpiece</p> <p>1-21 HZ</p> <p>Max 5.0 Bar</p>	<p>C15, DI15, D20-T Transmitter Heads</p> <p>Complete set of 9 transmitter heads</p> <p>V-Actor Handpiece</p> <p>Fascia Transmitter Heads</p> <p>Spinal Transmitter Heads</p>

Please note that optional accessories are highlighted in red text. Everything in black text is included with the device.



# MP 200

SIZE (WIDTH/HEIGHT/DEPTH)	KEY FEATURES	COMPATIBILITY
<p>46.6cm/18.7cm/45.4cm</p> <p>50lbs</p>	<p>Patented Handpiece Control Technology</p> <p>Falcon Handpiece</p> <p>1-21 HZ</p> <p>Max 5.0 Bar</p> <p>Vacu-Actor Cupping Technology and Materials Included</p> <p>Storz Tablet Included</p>	<p>C15, DI15, D20-T Transmitter Heads</p> <p><b>Complete set of 9 transmitter heads</b></p> <p>Storz Tablet</p> <p>Vacu-Actor Handpiece</p> <p><b>Fascia Transmitter Heads</b></p> <p><b>Spinal Transmitter Heads</b></p> <p><b>Storz Cart</b></p>

Please note that optional accessories are highlighted in red text. Everything in black text is included with the device.